



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



KANGAROO 924 enteral feeding pump (Ex-Display)

New (post-show),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
U.K. production.
Power supply: 230 V,
Emergency power supply: via 1 built-in new 7.2 V / 1.3 mAh battery,
Frequency: 50/60 Hz,
Enteral feeding pump
LP KANGAROO 924
constitutes a rotary peristaltic pump designed to control the flow rate of enteral feeding products,
Features:
Built-in rechargeable battery for portability,
Wide range of flow rates: 1-300 ml/hour, adjustable in 1 ml/hour increments,
Compact size and low weight for easy transportation,
Operation via touch panel and dial,
Automatic component check when the pump is turned on,
Large liquid crystal display,
Low occlusion pressure - the pump will not overcome back pressure of more than 12 psi (82.7 kPa),
Adjustable alarm volume,
Indicator lights for a functioning pump,
Memory of flow rate settings is retained for a period of 16 or 24 hours after the pump is turned off, depending on the state of connection to the AC source,
Audible and visual alarms:
The feed container is empty,
The feed tube becomes clogged,
Droplet sensors are blocked ,
The battery is discharged (the alarm goes off 15 minutes before the pump is turned off) ,
The device has been paused for more than 5 minutes,
Improperly loaded feeding set,
A self-detected problem requiring technical service has occurred ,
Pump motor failure occurred,
Ambient lighting error,
Included:
User manual in ENG, FR, GER, IT, ESP, SWE, DUT, PT, FI, DAN, GREK, CS, HUNG, RUS, PL, TUR, NOR (PDF),
Power cord,
Dimensions: 14 x 10.8 x 18.4 cm,
Weight: 1.8 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

