



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Cryotherapy apparatus KRIOPOL R (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Polish production, (2016),
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 500 W,
KRIOPOL R
(version 1) is a mobile and easy-to-use device that allows you to carry out an effective and efficient cryostimulation procedure. The working medium in it is liquid nitrogen. Thanks to it, we obtain a temperature of -1600 C at the nozzle outlet, already after about 30 seconds after switching on the device, which provides effective and rapid stimulus stimulation,
The intensity of the blowing is adjustable in steps (4 stages), depending on the size of the area to be cooled. The temperature of nitrogen vapor is -130oC (10 cm from the nozzle outlet),
The use of cryotherapy treatments is justified due to the beneficial effects of cold on the human body. Cryotherapy inhibits inflammatory processes by slowing down enzymatic reactions,
The beneficial effects of low temperatures are also manifested by a decrease in muscle tension, a decrease in the rate of conduction in nerve fibers, inhibition of nociceptors in the skin, blocking C-fiber and reducing the release of pain mediators,
Cryostimulation, causing a number of positive reactions in the body significantly increases the effectiveness of kinesitherapy performed immediately after cryotherapy,
Technical specifications:
Capacity: 30 l.,
Refrigerant: liquid nitrogen (LN2),
Gas jet temperature: -160°C (at the nozzle outlet) ,
Number of degrees of blowing intensity control: 4,
Liquid nitrogen consumption: 0.08 to 0.15 kg/min (continuous operation), depending on the set blowing intensity,
Number of treatments from one filling of YDS 30 tank: 40 to 60 (for average treatment time = 3 min),
Included:
Power cord,
Dimensions: 47 x 54 x 100 cm,
Weight: 34.5 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

