



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **OLYMPUS KV-5 endoscopic suction machine (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
U.K. production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 130 VA,  
Simple and intuitive operation,  
Suitable for all endoscopic procedures,  
Variable suction with vacuum pressure gauge: 0 to -1 bar,  
Maximum suction pressure: -85 kPa,  
Quiet pump operation,  
Includes:  
Power cord,  
User manual in English (PDF),  
Seeres jar with 1000 ml disposable cartridge,  
Hydrophobic microbial filter,  
Dimensions: 27 x 39 x 22 cm,  
Weight: 13 KG,  
It has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

