



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ERBE ERBOKRYO AE Cryosurgery Apparatus (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 40 W,
Application:
Ophthalmology,
Dermatology,
Otolaryngology,
Proctology,
Veterinary Surgery,
Plastic surgery,
Gynecology,
Thoracic surgery,
Specifications:
Temperature setting range: -20 °C to -70 °C,
Defrost temperature: From 8 °C to 15 °C,
Cryogenic gas:
Gas type: Nitrous oxide (N₂O) or carbon dioxide (CO₂),
Operating pressure: 40 - 60 bar,
Gas consumption "freezing": 36 g/min,
Amount of gas produced: 40 - 60 l/min,
Instrumentation placed on ERBE mobile table:
Mobile 4-wheel (2 braking wheels),
Pull-out drawer for accessories: 48 x 32 x 8 cm,
Built-in cylinder mount,
Includes:
ERBE 20716-087 Ø 3.0 mm probe,
Dedicated ERBE mobile table,
Hose for supply cylinder,
Power cord,
Dimensions:
Appliances: 47 x 35 x 10 cm,
Overall: 64 x 50 x 90 cm,
Weight:
Devices: 10 KG,
Total: 25 KG,
It has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

