



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **VWR INCU-LINE IL23 Incubator (Ex-Display)**

New (post-show),

Real photos of the product,

Slovenian production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 160 VA,

Digital incubator

INCU-LINE IL23

offers an economical and space-saving solution for microbiology or hematology applications,

The temperature can be freely adjusted in 0.1 ° C increments up to a maximum of 70 ° C,

Heating elements are built into the base and walls of the enclosure to ensure temperature stability and uniformity,

The housing and door frame are made entirely of epoxy-coated mild steel,

Technical specifications:

Ambient From 5°C to 70°C,

Total capacity: 23 l,

Air circulation: natural convection,

Included:

Power cord,

2 pcs stainless steel shelf 28.5 x 28 cm,

Dimensions:

Internal: 29 x 29 x 27,

External: 34 x 35 x 40 cm,

Weight: 12.8 KG,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

