



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Psychological apparatus for measuring drivers' reaction time ALFA ELECTRONICS (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Polish production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Application:  
Drivers,  
Candidates for drivers,  
Examiners,  
Candidates for instructors and examiners,  
Operators (forklifts, cranes, overhead cranes, excavators, etc.),  
Athletes,  
In other scientific fields (such as ergonomics),  
Presentation of results:  
Average reaction times,  
Minimum time,  
Maximum time,  
Range,  
Number of wrong reactions,  
Number of missed responses,  
Number of reactions to the wrong stimulus,  
Number of reactions with the wrong button,  
Number of responses between stimuli,  
Advantages:  
Very simple operation,  
Standardized test programs permanently stored in the device,  
Aesthetic easy-to-clean housing,  
Small size and weight for portability and exit examinations,  
Cooperation with computer and ODOIA program,  
Menu in language: Polish,  
Included:  
Power cord,  
Dimensions: 60 x 45 x 96 cm,  
Weight:  
Has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

