



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Audiometer ESSILOR AUDIOSCAN OTO (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real product photos,

French production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Technical specifications:

Built-in 5.9" LCD display,

Frequency range: 125 - 16,000 Hz,

Smooth volume control: From -10 to 100 dB,

Programmed modes:

Audioscan test,

BEKESY test,

ISCAN test,

H & W test,

AC MANUAL test,

BC MANUAL test,

Strategies,

Built-in thermal printer with a paper width of 11.5 cm,

Possibility to select the menu in the following languages: FR, ENG, IT, ESP,

Included:

Power supply with cable power supply,

Headphones,

Patient button,

Dimensions: 22 x 35 x 12 cm,

Weight: 4 KG,

Cat 02,

It has a current inspection and is ready to work,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,

Financing options (Poland only): Installments, Leasing, Loan,

