



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



KARL STORZ TELECAM PAL camera processor 20210020 (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Manufactured by U.S.A.,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 15 W,
Specifications:
White balance function,
Automatic exposure control mode,
Manual exposure control mode,
Features the following connectors:
2x S-VIDEO,
GENLOCK IN,
KEYBOARD IN,
OUTPUTS:
GREEN,
RED,
BLUE,
SYNC,
COMP VIDEO

Included:
Power cord,
Dimensions: 30 x 32 x 8.5 cm,
Weight: 3 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

