



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Portable low-temperature freezer STIRLING ULTRACOLD ULT-25NEU (Ex-Display)

New (post-show),

Real photos of the product,

Manufactured by U.S.A., (2019),

Power supply: 230 V / 12 V DC,

Frequency: 50/60 Hz,

Power: 280 W,

STIRLING ULTRACOLD ULT-25NE

is a portable low-temperature freezer that offers compact, reliable, energy-efficient storage and transport of biological materials at temperatures from -20°C to -86°C. Its small size and light weight make this freezer model ideal for anyone transporting materials that require special storage conditions,

Specifications:

Capacity: 25 liters,

Easy to transport,

Small size,

Lockable compartment,

Cooling temperature: from -20°C to -86°C,

Uses 100% natural refrigerants,

Includes:

English instruction manual (PDF),

Power cord,

12 V power cord (cigarette lighter socket),

Key (1 pc.),

Basket (2 pcs.),

Dimensions: 70 x 36 x 46 cm,

Weight: 21 KG,

It has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

