



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



VIASYS HEALTHCARE MICROLAB Spirometer (Reconditioned)

Reconditioned (used),
Technical condition: good,
Visual condition: good,
Real photos of the product,
U.K. production,
Power supply: 230 V,
Emergency power supply: via new built-in 8.4 V battery,
Frequency: 50/60 Hz,
The MicroLab VIASYS portable spirometer is equipped with a high-quality touch screen and a fast and quiet printer,
Using easy-to-read icons (corresponding to available functions) and a stylus or optional mouse, operation of the spirometer is extremely simple,
Unique context-sensitive help available at the press of a single key explains the operation of all spirometer functions,
Technical specifications:
Software: V3.02,
Memory: 2000 tests (flow/volume loop tests, volume/time curves),
Data entry: first name / last name / ID number / age / weight / height / gender,
Display: 1/4 VGA color LCD,
Battery operation for 30 hours on a full charge,
Features:
Reliable turbine - accurate results especially at low flows,
Test validation,
Comparison of results with due values,
Supplied in an aesthetic case with accessories,
Lung age and interpretation according to Enright, ATS or BTS,
Memory for 2000 tests,
41 parameters measured, Inspiratory and expiratory spirometry, intensity and free spirometry,
Calculation of body mass index (BMI) and dyspnea scale,
Built-in high-resolution printer allows printing of flow/volume and volume/time curves,
Tests after administration of bronchodilator and/or steroid,
No calibration required,
Selectable menus in CN, DE, EN, ES, FR, IT, DK, PL, PT, SE, RU,
The device has traces of use that do not interfere with the comfort of use,
Included:
User manual in English (PDF),
Power cord,
2 pcs. - turbines,
Head,
Power supply,
Mouthpiece - 8 pcs,
Dedicated carrying case,
Dimensions: 12 x 25.5 x 2.5 cm,
Total weight: 1 kg,
It has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

