



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



INTERACOUSTIC AD229e diagnostic audiometer (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Power supply: 230 V,
Frequency: 50/60 Hz,
High-quality diagnostic audiometer type
AD229e
works well for both stationary and portable testing,
Technical specifications:
Level intensity: AC: From -10 to 120 dB HL - gradual adjustment in 1 dB / 5dB increments,
BC: From -10 to 80 dB - gradual adjustment in 1 / 5 dB increments,
Frequency range: 125-8000 Hz,
Masking: Narrowband(or white noise) for sound and speech presentation,
Communication: built-in microphone "talk forward" volume control 0-110dB SPL,
Tone presentation:
Manual,
Reverse "Reverse",
Single pulse,
Multi pulse 250-5000 msec.
On / Off,
Frequency selection: 125Hz, 250Hz, 750Hz, 1500Hz / 8kHz,
Additional tests:
ABLB. Stenger,
Stenger Speech,
HW automatic test,
SISI,
Langenbeck (tone in noise). Békésy,
2 channel speech test,
2 channel Master Hearing Aid,
Cooperation with software:
Effetha Audio,
Interacoustics OtoAccess™,
NOAH,
Control panel keyboard in language: English,
Included:
User manual in English (PDF),
Power cord,
TELEPHONICS in-ear headphones,
Patient button,
Dimensions: 36 x 26 x 10 cm,
Weight: 2.5 KG,
It has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

