



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NOUVAG AG DISPENSER DP 30 tumescent anesthesia pump (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Swiss production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 40 VA,
Dispenser DP 30
is a universal pump used for tumescent anesthesia in endovascular, liposuction and other procedures,
Intuitive to use,
Equipped with a metal stand for hanging a bottle of tumescent fluid,
Infinitely precise flow rate control with potentiometer from 0 - 100% ,
Foot switch allows you to turn the flow on and off and control the discharge speed without interrupting the operation,
Specifications:
Flow rate: 1 - 15 L/h,
Maximum pressure: 2 bar,
Includes:
38 cm high tripod for hanging the tumescent fluid bottle,
NOUVAG Model VARIO foot attachment,
Multi-language instruction manual (PDF),
Power cord,
Dimensions: 25.5 x 23 x 48 cm,
Weight: 3.5 KG,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

