



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MIDMARK gynecological chair 604-004 (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Manufactured by U.S.A,  
MIDMARK RITTER 304  
is a gynecological examination chair, designed for both offices and universities,  
Upholstery with medical certificate in red color,  
Smooth manual adjustment:  
Backrest with spring release handle, located on the left and right side of the table in the range: 0° - 70°,  
Seat: 0 - 5°,  
Fixed level recliner height: 86 cm,  
Features:  
Extendable, foldable and side-adjustable stirrups,  
Removable waste bowl, made of plastic with dimensions: 23.5 x 29 x 7 cm,  
A retractable platform for easy access to the chair,  
Armchair extension with dimensions: 46 x 50 x 5.5 cm,  
Two front drawers for additional accessories in size: 43 x 48 x 10 cm,  
Three side drawers with dimensions: 54 x 37.5 x 10 cm,  
Dimensions:  
Backrests: 68 x 88 cm,  
Seats: 72 x 50 cm,  
Length in reclining position: 144 cm,  
Overall: 74 x 140 x 83/165 cm,  
Maximum load capacity: 227 KG,  
Total weight: 100 KG,  
Has a current inspection and is ready for use,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

