



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SCHMITZ MEDI-MATIC gynecological chair 115.320 (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
New medical approved upholstery in the color of your choice from our color chart,
Schmitz Medi - Matic
Used for gynecological and urological examinations and procedures,
Technical specifications:
Electromechanical adjustment with foot control:
Manual adjustment of 2-segment back rest in the range: 0° / 42°,
Electric adjustment:
Height: 84 - 104 cm,
Slant adjustment: +6° / 0° / -30°,
Accessories:
Double-sided rails for attaching additional accessories,
Pull-out patient step,
Pull-out waste bowl with dimensions: 32 x 30 x 6 cm,
Built-in 2 outlets for powering additional equipment,
Built-in storage compartment for paper roll,
Double-sided leg support,
Removable leg rest,
Includes:
Power cord,
Wired foot controller,
2 pcs. (leg rest),
Leg rest,
Dimensions:
Headrest: 60 x 22 x 10 cm,
2-segment backrest: 60 x 64 x 10 cm,
Seat: 60 x 32 x 10 cm,
Leg rest segment: 62 x 67 x 10 cm,
Total: 70 x 105 x 150 cm,
Maximum load capacity: 135 KG,
Weight: 173 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

