



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MEMMERT BE 600 laboratory incubator (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 1600 W,  
Technical specifications:  
Infinitely adjustable temperature: from 20° to 70°C,  
Infinitely adjustable time: from 1 minute to 999 hours or continuous operation mode,  
Chamber interior made of stainless steel,  
Chamber capacity: 256 l,  
Chamber equipped with two sieve trays with possible arrangement,  
Over-temperature alarm,  
Built-in digital temperature display,  
Adjustable fresh air intake,  
Additional internal glass door,  
Built-in RS 232 serial interface,  
Included:  
Power cord,  
Dimensions:  
Overall: 95 x 68.5 x 91.5 cm,  
Internal: 80 x 50 x 64 cm,  
Weight: 87 KG,  
Has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for a period of 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

