



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ACUTRONIC MONSOON III nozzle ventilation ventilator (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Swiss manufacture, (2018),
Power supply: 230 V,
Emergency power supply: through new built-in batteries,
Frequency: 50/60 Hz,
Power: 250 W,
Large 9-inch touchscreen displaying patient parameters, graphical waveform, pressure levels and limits,
Greater safety and reliability with microprocessor control,
Adjustable bypass flow allows pre-oxygenation of patients before surgery and auxiliary ventilation,
Features:
Large 9-inch detachable touchscreen display showing patient parameters, graphic waveform, pressure levels and limits,
Displays delivered tidal volume and minute volume,
Adjustable bypass flow for pre-oxygenation of patients before surgery and auxiliary ventilation,
Greater safety and reliability with microprocessor control,
PIP and PP alarm limit/high pressure cutoff,
MAP, EEP, PIP, PP pressure monitoring displayed on screen for easy viewing,
Displayed and measured settings:
Frequency,
Inspiration Time,
Pause Pressure,
Drive Pressure,
FiO2,
TV / MV,
PIP / MAP / EEP,
Humidification,
Bypass flow,
Alarms:
Pause Pressure (PP),
Peak pressure (PIP),
OFF. - 101 cmH2O,
1 - 101 cmH2O,
Parameter settings:
Frequency: 12 - 150 CPM,
Inspiration time: 20 - 60% ,
Fio2: 21 - 100% oxygen,
Drive pressure: 5 - 40 PSI,
Auxiliary gas flow rate: 0 - 70 lpm,
Humidification: 0 - 100% ,
Instrumentation placed on PRO-CART mobile cart,
Selectable menu in GER, ENG, FR, IT, SWE, PL, HUN, RUS, CS,
Included:
Power cord,
Gas supply air lines,
Can be retrofitted with consumables at additional cost,
Dimensions:
Respirator: 42 x 39 x 19 cm,
Total: 46 x 70 x 105 cm,
Total weight: 40 KG,
Has up-to-date inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

