



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



KARL STORZ endoscope cart 29003BH (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
KARL STORZ 29003BH
Is a stable as well as sturdy endoscope cart with freedom of movement,
It has:
Built-in transformer,
Top and middle table top with dimensions: 43 x 48 cm,
Lower top with dimensions: 50 x 48 cm,
Pull-out drawer with dimensions: 39 x 38 cm,
Tripod with free height adjustment,
3 built-in power cords for additional accessories,
Rubber bumpers,
Mobile cart on 4 wheels (2 braking),
Includes:
Power cord,
3 power cords additional accessories,
2x key for pull-out drawer,
Pull-out tripod,
Inbus key,
Overall dimensions: 54 x 66 x 138 cm,
Maximum load capacity:
Top shelf and middle shelf: 50 KG,
Lower shelf: 225 KG,
Total weight: 38 KG,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

