



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Infusion (syringe) pump B. BRAUN PERFUSOR COMPACT (Cat 10) (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Power supply: 230 V,
Emergency power supply: via 4 new built-in 4 x 1.2 V AA batteries,
Frequency: 50/60 Hz,
Power: 10 VA,
Perfusor Compact
is a portable syringe pump, prepared for intravenous fluid therapy, nutrition,
Simple and quick exchange of 20 and 50 ml syringes from different manufacturers and their reliable attachment with automatic volume recognition,
Easy selection of functions, including changing the speed without interrupting the infusion,
Sensitive, adjustable pressure sensor provides a high level of safety with different dosing modes, when the system pressure increases critically, the infusion speed automatically decreases,
Device parameters:
Operating time:
Approximately 10 hours at a flow rate of 99.9 ml/h,
Approximately 80 hours at an average flow rate of 10 ml/h,
Choice of syringes: 20 or 50 ml,
Delivery rate: from 0.1 to 99.9 ml/h (in 0.1 ml increments),
Bolus rate: 800 ml/h,
Infusion volume: from 0.1 to 999.9 ml (in 0.1 ml increments),
Features:
Foiled control panel,
Liquid crystal display,
Menu in Polish,
Audible and visual alarms,
Multifunctional MFC connection,
Comfortable carrying handle,
Possibility to assemble the kit,
Supported syringes:
20 ml:
Monoject EU (1100-620036),
Monoject USA (8881-520657),
B-D Pipak,
Terumo (SS 20L),
B|Braun Omnifix LL,
Original Perfusor,
50/60 ml:
Monoject EU (1100650090),
Monoject USA (8881-560125),
Terumo (BS5OLG),
B-D Pipak,
50 ml:
Original Perfusor,
B|Braun Proinject,
B|Braun OMNIFIX LL,
B|Braun Euroject,
60 ml:
Terumo USA (SS 60LL),
Includes:
User manual in Polish and English (PDF),
Power cord,
Dimensions: 24.5 x 14.5 x 10.5 cm,
Weight: 1.4 kg,
Cat 10,
Has an up-to-date inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

