



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Wheelchair WZSO (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Polish production,
WZSO wheelchair, thanks to its low weight, the use of crossbar design, it can be easily folded and transported in the car,
Technical specifications:
Foldable lightweight steel frame with cross-brace construction,
Ergonomic handles for the caregiver,
Footrests equipped with additional calf support belt,
Brake for the patient,
Maintenance-free front wheels,
Rear wheels pumped with a diameter of 24" inches,
Dimensions:
Seats: 42 x 41 cm,
Backrests: 42 x 33 cm,
Overall: 70 x 90 x 90 cm,
Weight:
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

