



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Patient lift verticalizer VERMEIREN ALBATROS (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Belgian production, (2021),  
Power supply: 230 V,  
Emergency power supply: through 2 built-in 12 V / 2.9 Ah batteries,  
Frequency: 50/60 Hz,  
Albatros Vermeiren  
is an uprighting device. It makes it possible to verticalize a patient from a sitting to an upright position. The device is also used to transfer people from a wheelchair or bed to the toilet in an upright position,  
The Albatros Lift supports caregivers and medical personnel in everyday situations and effectively relieves the backbone during transfers of disabled people thus reducing the risk of musculoskeletal injuries. With this device, it is possible for only one caregiver to perform such activities,  
Product Features:  
Four points of support : feet, knees, hips, chest,  
Entire control and drive system based on LINAK technology,  
Design allows easy and quick folding for transportation without using any tools,  
Models equipped with locking casters,  
Suitable foot platform equipped with non-slip inserts and proper knee support with calf protection dog perfectly stabilize the legs. If more foot stabilization is required, special foot troughs with stabilizing straps can be used,  
The running frame is movable and allows smooth spacing adjustment for maximum stability,  
Specifications:  
Lowest arm position: 104 cm,  
Highest arm position: 173 cm,  
Arm lifting range: 70 cm,  
Built-in central brake,  
Electric control with wired remote control,  
Includes:  
User manual in Polish (PDF),  
Power cord,  
2 pieces of patient carriers,  
Wired control remote control,  
Overall dimensions: 70 x 110 x 170 cm,  
Maximum lifting capacity: 150 KG,  
Weight: 45 KG,  
It has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

