



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



BOBBY VERMEIREN EVO wheelchair (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Belgian production, (2020),
Wheelchair
Bobby
is a professional transportation equipment, dedicated to people who have problems with free movement, It is designed for indoor and outdoor use,
Wheelchair
Bobby
is an extremely lightweight and easy to fold product. It is ideal for use in hospitals, airports, and transportation from point A to point B. What's more, the transport equipment is convenient and quick to fold, we can easily store it in very small spaces and a small car trunk,
Features:
Quick folding design option with convenient folding backrest,
Equipped with brakes for the companion,
Features soft and comfortable armrests,
Shallow footrest with adjustable angle,
Additional belt to prevent the patient from sliding off the seat,
Practical handles on both sides for carrying the wheelchair,
Modern and elegant design,
Technical specifications:
Minimum turning radius: 150 cm,
Seat inclination: 5°,
Hand brake on both sides,
Rear wheel diameter: 12.5",
Dimensions:
Seat: 48 x 40 cm,
Backrest: 40 x 42 cm,
Overall: 60 x 104 x 92 cm,
Weight: 13 KG,
Has an up-to-date inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

