



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Operating table STERIS CMAX-T2L (Ex-Display)

New (exhibition),
Real product photos,
French production,
Power supply: 230 V,
Emergency power supply: via 2 new built-in batteries 12 V / 24 Ah,
Frequency: 50/60 Hz,
Power: 440 VA,
STERIS CMAX-T2L mobile operating table with electro-hydraulic drive,
Emergency control panel built-in in the table column,
The entire 5-segment table top is transparent to X-rays,
Mounted side rails on both sides,
Built-in side rails for mounting additional equipment,
Medically certified upholstery in black (possibility of replacing the upholstery with a new one for an additional fee),
Adjustable in many planes using a wireless remote control or via a built-in panel on the column:
UP/DOWN - operating surface height (up/down): 68 - 114 cm,
TREND/REV - changing the patient's lying direction relative to the table column: Trendelenburg position (30°) and Reverse Trendelenburg (44°),
TILT L/TILT R - lateral tilt in the range of 25°,
LEG UP/LEG DOWN - leg segment adjustment: up/down in the range of: -90° / 0° / 90°,
LEVEL - zeroing: returning to the horizontal position of the bed,
Longitudinal shift - front / back: 22 cm,
FLEX / REFLEX - possibility of bending the table top in both directions: roof position / V position (Beach chair),
Built-in electric wheel lock,
Manual adjustment:
Headrest tilt: 45° / 0° / -40°,
2-segment leg rest:
Upper part - thigh rest: -90° / 0° / 90°,
Lower part - calf rest: -90° / 0° / 90°,
Free adjustment of the leg rest opening in the range of: 45°,
Possibility of extending the calf rest in the range of: 33 cm,
Included:
Power cable,
Wireless (Bluetooth) handheld remote control,
Dedicated cable for charging the handheld remote control,
Dimensions:
Headrest: 48 x 28 x 6 cm,
Bed segment: 84 x 52.5 x 6 cm,
Thigh segment: 24 x 30 x 6 cm,
Calf segment: 24 x 44 x 6 cm,
Total length:
Maximum load capacity: 270 KG,
Weight: 270 KG,
It has a current inspection and is ready to work,
Issued Technical Passport (Service Report) valid for the period 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international market,
Financing options (Poland only): Installments, Leasing, Loan,

