



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Digital flat scale SECA 861 (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- Power supply: 230V,
- Emergency power supply via new built-in 9 V batteries,
- Frequency: 50/60 Hz,
- The Seca 861 digital flat scale is a practical scale, equipped with a clear LCD screen that turns on and off automatically, and its batteries ensure a long battery life of at least 16,000 weighings,
- Stability and durability combined with high autonomy make the SECA 861 electronic digital floor scale the ideal choice for everyday, intensive use in clinics and healthcare facilities with high traffic,
- **The SECA 861** has a raised rubber surface that ensures stability and is easy to clean,
- **Specification:**
 - Graduation: 50 g <150 kg> 100 g,
 - Scale type: Flat scale,
 - Scale Type: Digital Function,
 - Automatic shutdown,
- **Included in the set:**
 - Power cord,
- Dimensions:
 - Platforms: 27 x 27 cm,
- Maximum load capacity: 200 KG,
- Weight: 8 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for domestic market (Poland),**
 - **3 months to international market,**
 - **Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,**
- Financing options (Poland only): Installments, Leasing, Loan,

