



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ENRAF NONIUS CURAPULS 419 short wave diathermy (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Dutch production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 400 W,
Technical specifications:
Frequency of emitting waves: $F = 27.12$ MHz,
Foiled control panel,
Operating (wave) mode: continuous or pulsed,
Infinite power adjustment: 0 - 10 (dial),
Timed operation switch (electronic) in the range of 0 - 30 min (increments of 0.5 min),
Gradual adjustment of Hz pulse frequency: 15 / 20 / 26 / 35 / 46 / 62 / 82 / 110 / 150 / 200 - dial,
Applicator placed on separate adjustable arms in multiple planes,
Included:
User manual in ENG, FR, NED, GER, ES (PDF),
Power cord,
Applicators (induction electrode): diameter (\varnothing) 14 cm,
Overall dimensions: 80 x 44 x 95/153 cm (height with arm dissolved),
Weight: 68 KG,
It has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

