



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **ENRAF NONIUS CURAPULS 419 short wave diathermy (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Dutch production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 400 W,  
Technical specifications:  
Frequency of emitting waves:  $F = 27.12$  MHz,  
Foiled control panel,  
Operating (wave) mode: continuous or pulsed,  
Infinite power adjustment: 0 - 10 (dial),  
Timed operation switch (electronic) in the range of 0 - 30 min (increments of 0.5 min),  
Gradual adjustment of Hz pulse frequency: 15 / 20 / 26 / 35 / 46 / 62 / 82 / 110 / 150 / 200 - dial,  
Applicator placed on separate adjustable arms in multiple planes,  
Included:  
User manual in ENG, FR, NED, GER, ES (PDF),  
Power cord,  
Applicators (induction electrode): diameter ( $\varnothing$ ) 14 cm,  
Overall dimensions: 80 x 44 x 95/153 cm (height with arm dissolved),  
Weight: 68 KG,  
It has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for a period of 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

