



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **LEICA M300 dental microscope (Reconditioned)**

Product temporarily unavailable  
Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Halogen bulb: 100 W,  
Specifications:  
Magnification: 6.4x / 10x / 16x / 25x / 40x,  
Eyepiece: 10x / 21,  
Freely adjustable eyepiece spacing,  
Light intensity adjustment with a knob,  
Height adjustment of the arm in the range: 55 - 158 cm,  
Hand mount for free adjustment of microscope rotation,  
Microscope equipped with light source BOB Fiber Optic Technique,  
Microscope placed on a mobile 4-wheel base (2 braking wheels),  
Included:  
Power cord,  
Dimensions:  
Base diagonal: 60 cm,  
Total height of the top of the arm: 218 cm,  
Weight: 40 KG,  
Has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,  
The warranty period does not include bulbs and other light sources installed in the device,  
Possibility of financing (Poland only): Installments, Leasing, Loan,

