



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Medical assistant (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- Mobile medical assistant with pull-out drawers,
- Pull-out top shelf dimensions: 38 x 41 cm,
- 3 pull-out drawers with dimensions: 35 x 36 x 9 cm,
- Mobile assistant on wheels,
- Overall dimensions: 45 x 48 x 78/90 cm,
- Weight: 10 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
  - **6 months for domestic market (Poland),**
  - **3 months to international market,**
  - **Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,**
- Financing options (Poland only): Installments, Leasing, Loan,

