



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STRYKER Arthroscopic Kit (SV-2, L9000, 1088i HD, FLUID SAFE) with STRYKER HD Camera Head (1088-210-122) (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
The set includes:
STRYKER SV-2 240-030-920 medical monitor:
Made in the U.S.A.,
Power supply: 230 V,
Frequency: 50/60 Hz,
Specification:
Screen diagonal: 19" inches,
Matrix: LCD,
Resolution: 1280 x 1024 px.,
Response time: 8 ms,
Connectors: VGA / DVI / C-Video / S-Video / RS-232,
STRYKER L9000 light source:
Made in the U.S.A.,
Power supply: 230 V,
Frequency: 50/60 Hz,
Stryker light source The L9000 is an advanced lighting system designed to enhance visibility during minimally invasive surgical procedures.

Using LED technology, it delivers bright, clear light via fiber optics to provide clear visualization of surgical sites for arthroscopy, laparoscopy, and general endoscopy applications.

Features:
Bright, clear LED illumination: Provides consistent, high-quality light for clear visualization during surgery.

Electronic Scope Sensing Technology (ESST): Automatically switches the device to standby mode when the cable is not connected to the scope, increasing safety.
Multiple Procedure Compatibility: Ideal for arthroscopy, laparoscopy, and general endoscopy across a variety of surgical specialties.

Universal Connectivity: Compatible with all Stryker light cables and can be adapted to other endoscopes with the appropriate adapters.
Efficient LED technology: Reduces heat generation while providing long-lasting, reliable performance during procedures,

Specification:
LED technology,
Free 10-step light intensity adjustment,
3.5" touch screen,
STRYKER 1088i HD image processor:
Made in the U.S.A.,
Power supply: 230 V,
Frequency: 50/60 Hz,
STRYKER 1088i HD camera processor is a device used to record high quality video images in endoscopic and surgical applications,
The front control panel has built-in functions:
White balance operation,
4-step light intensity adjustment,
Increasing or decreasing image sharpness,
Possibility to select surgical specialization,
Image processor equipped with STRYKER HD 24 mm camera head (1088-210-122),
STRYKER FLUID SAFE arthroscopic pump:
Made in the U.S.A. U.S.A.,
Power supply: 230 V,
Frequency: 50/60 Hz,

Stryker FluidSafe is a system used in medicine, especially in surgery, to safely collect and dispose of body fluids (e.g. blood, tissue fluid, rinsing fluid) during surgical procedures,

Features:
Fluid suction – used to effectively remove fluids from the surgical field, maintaining cleanliness and visibility for the surgeon,
Fluid filtration and collection – the system can filter and collect fluids to enable their later disposal in accordance with medical regulations,
Minimizing the risk of infection – prevents the spread of body fluids, which reduces the risk of cross-infection,

Specification:
Maximum output pressure: 200 mm Hg,
Pressure measurement accuracy: ±1.5 %,
Maximum flow: 500 ml/min.,

SONY RDR-HX720 DVD burner:
Specification:
Disk capacity: 160 GB,
Media recording: CD / CD-R/W / DVD+R/W / DVD-R/W,
Video file reading: DVD-Video / DivX,
Audio file reading: MP3,
Image file reading: JPEG,
Video processing: 12-bit 108Mhz,
Audio processing: 24-bit 192khz,
Equipment placed on the mobile ALVO MEDICAL arthroscopic tower:
Double-sided handles for easy movement,
Built-in 4 wheels (2 braking),
Shelves with dimensions: 55 x 60 cm,
Sliding drawer for accessories 49 x 55 x 10 cm,
Included:
Power cables,
Set of connecting cables,
Remote control for DVD recorder,
Overall dimensions: 64 x 60 x 196 cm,
Total weight with equipment: 130 KG,
It has a current inspection and is ready to work,
Issued Technical Passport (Service Report) valid for the period 12 months,
Warranty:
6 months for domestic market (Poland),
3 months for international market,
Possibility to extend warranty for an additional fee to 12 months for domestic and international market,
Financing options (Poland only): Installments, Leasing, Loan,
Reconditioned (used)

