



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



VITEA CARE walker (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

Polish production, (2020),

Walker

VITEA CARE

is perfect for any situation perfectly facilitating verticalization and provides support for people with reduced digestibility of movement,

It has a soft fabric seat that allows you to rest comfortably,

The four-wheeled structure is entirely made of lightweight aluminum, which provides greater comfort,

The two-sided handles have anatomical shape and are made of soft plastic, which increases the comfort of its use,

The walker has a fabric container for essentials,

Dimensions:

Height in range: 78 x 93 cm,

Seat height, 51 - 57 cm,

Seat width: 35 cm,

Maximum load capacity: 120 KG,

Weight: 7 KG,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

