



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MAQUET RADIUS gynecological chair 1557.01AJ (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 290 VA,  
MAQUET RADIUS  
is used for gynecological and urological examinations and procedures,  
The upholstery will be replaced with a new one in a color chosen by the customer from our color chart,  
Electromechanical adjustment by foot attachment and wired handheld remote control:  
Seat height in the range of 60 - 90 cm,  
Tilt angle:  
Back segment in the range of (- 9°) / 0° / 57°,  
Seat in the range of 18°,  
When adjusting the angle of the back segment, the seat angle is automatically adjusted,  
3 preset buttons with function memory,  
STOP function - button with activation of immediate interruption of the function being performed,  
The chair has:  
Storage compartment for a 38 cm wide roll of backing paper,  
A cutout in the upholstery of the back segment for the backing paper,  
Retractable waste shelf measuring 18 x 25 cm,  
Reversible footrests,  
Removable protective film for the footrests,  
Ability to unfold the chair to the reclining position,  
Includes:  
Foot attachment,  
Wired remote control,  
Waste bowl,  
Dimensions:  
Back segment: 91 x 61 cm,  
Seat: 61 x 31 cm,  
Length in extended position: 175 cm,  
Overall: 158 x 75 x 138/168 cm,  
Maximum load: 180 kg,  
Weight: 108 kg,  
Has a current inspection and is ready for use,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

