



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ARTROMOT TRAPO K1 mobile stand for CPM rail. (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

German manufacture,

Mobile cart

ARTROMOT TRAPO K1

for CPM rail with built-in adjustment of bilateral tilt to facilitate the performance of rehabilitation,

The wheelchair has 4 wheels with built-in brakes and has an adapted frame for performing rehabilitation on a patient in the supine position,

It has a single-stage 90° tilt on both sides,

4 material fasteners along with anti-slip on both sides allowing the rail to be firmly fixed,

Dimensions:

Rail top: 40 x 31 cm,

Overall: 72 x 70 x 77 cm,

Total weight: 10 KG,

Has an up-to-date inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

