



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



EPPENDORF CENTRIFUGE 5415D Laboratory Centrifuge (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 180 W,
Digital benchtop centrifuge
Eppendorf 5415D
is extremely easy to operate thanks to its knob control, and works practically like a ready-to-use machine,
Its small size and brushless motor make it maintenance-free. The 5415D's maximum speed is 13,200 rpm, and the digital timer allows it to run for up to 30 minutes,

Features:

Easy to use with dial control, ideal for routine tasks,
Ready-to-use, ready-to-use solution for convenience,
Maintenance-free thanks to brushless motor,
Maximum speed of 13,200 rpm for efficient spinning,
Digital timer with up to 30 minutes running time,
Space-saving, suitable for tabletop use,

Specifications:

Maximum speed: 13,200 rpm,
Maximum RCF: 16,110 grams,
Maximum capacity: 24 x 1.5/2.0 ml,
Timer: up to 30 minutes,
Acceleration time: less 13 seconds,
Deceleration time: less 11 seconds,
Centrifuge equipped with EPPENDORF F-45-24-11 rotor with dedicated lid,

Includes:

Instruction manual in English (PDF),
Power cord,
EPPENDORF F-45-24-11 rotor,
Dimensions: 23 x 31 x 23 cm,
Weight: 21 KG,

It has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

