



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



VITEA CARE EASY TRANSFER transport and toilet cart (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

Polish production, (2020),

Vitea Care Easy Transfer

is a professional assistive wheelchair that facilitates/replaces the use of the toilet and shower,

Easy Transfer

is great for transferring a patient from bedroom to bathroom for toileting, and can also serve as a typical indoor transportation equipment for disabled people to a limited extent,

Features:

White painted steel frame,

Height adjustment of the wheelchair,

Tube handle for mounting an IV, for example,

Seat height allows overrun on the toilet seat,

Closing the backrest with a strap with buckle and additionally with a screw,

Metal footrests, tiltable,

All solid wheels with brakes,

Dimensions:

Overall: 58 x 78 x 87.5 / 107.5 cm,

Seat height: 43 - 63 cm,

Armrest height: 68.5 x 88.5 cm,

Seat: 41 x 41 cm,

Maximum load capacity: 150 KG,

Weight: 19 KG,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,



