



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



2L KORGIEL W8 O2 + MT02 oxygen therapy kit (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

Polish production,

The reducer is used to reduce the pressure of gas taken from the cylinder to the level required by the user. The output pressure should be 4 or 8 bar depending on the consumer devices connected to the output,

Cylinder regulators work using a diaphragm, which gives very precise pressure regulation,

Oxygen dispenser with humidifier is used to dispense oxygen to the patient's airway and to humidify oxygen,

Technical specifications:

It is made of chrome-plated brass,

Built-in pressure gauge,

Pressure gauge indication range: 0-31.5 MPa,

Working pressure: 0.5 MPa,

Flow rate: 0 ÷ 17 l / min,

Humidification: 0 ÷ 60 % ,

Equipped with a special connection for oxygen cylinder,

Built-in rotametric dispenser,

Included:

Humidifier bottle:

Capacity: 200 ml,

Sterilizable (121°, time up to 20 minutes),

Thread: 9/16,

Gasket (oring),

UBT-01 bottle holder:

Adapted for oxygen cylinders with a diameter of 14 cm,

Cylinder attachment with buckled straps,

Mounting dimension: 35 x 54 cm,

Cylinder designed for medical oxygen with a capacity of 2 liters,

Reducer dimensions: 11.5 x 12 x 10 cm,

Dimensions of the set: 11.5 x 33 x 13 cm,

Weight: 1.1 KG,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

