



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Water Distiller MH943T-S (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 580 W,  
Water distiller of 4L capacity with stainless steel interior,  
Operation based on forced-air cooler,  
Application:  
For autoclave,  
Dental units,  
Air humidifiers,  
In dental offices,  
Laboratories,  
Specifications:  
Water tank capacity 4 liters,  
Capacity- approximately 1 liter per hour,  
Included:  
Power cord,  
Dimensions: 20 x 36 cm,  
Weight: 3 KG,  
Has a current inspection and is ready to work,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

