



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Water Distiller MH943T-S (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 580 W,
Water distiller of 4L capacity with stainless steel interior,
Operation based on forced-air cooler,
Application:
For autoclave,
Dental units,
Air humidifiers,
In dental offices,
Laboratories,
Specifications:
Water tank capacity 4 liters,
Capacity- approximately 1 liter per hour,
Included:
Power cord,
Dimensions: 20 x 36 cm,
Weight: 3 KG,
Has a current inspection and is ready to work,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

