



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LAMIDEY NOURY DUO MC2 (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
French production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 510 VA,
Surgical diathermy
LAMIDEY NOURY DUO MC2
is designed to perform outpatient and surgical procedures,
It has the following functions:
Monopolar cutting: 200 W. max,
Modes:
PURE,
HEMO 1,
HEMO 2,
Coagulation: 100 W max,
Modes
:
FORCEE,
DOUCE,
Infinitely adjustable power,
Digital displays of set power,
Foiled touch control panel,
Included:
Power cord,
Wired foot attachment,
Factory new cable with ERBE clip for neutral electrodes - (
No electrodes included, can be retrofitted at additional cost),
Factory new monopolar electrode cable with ERGO PEN BOWA JackKNIFE holder - 2 buttons, 2.4 mm shaft
(no electrodes included, can be retrofitted at additional cost)
,
Factory new bipolar cable, BOWA forceps
(can be retrofitted with bipolar tweezers at additional cost)
,
Dimensions: 42 x 42 x 18 cm,
Weight: 15 KG,
It has a current inspection and is ready to work,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

