



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Dental amalgam mixer amalgamator SYG-200 (Ex-Display)**

New (post-show),  
Real photos of the product,  
Manufactured by P.R.C., (2014),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 220 VA,  
SEPTODONT SYG-200 is a device aimed at obtaining an optimal mixture of biodentine, amalgams, cements and glasonomers.  
Specifications:  
Allows control to obtain a consistent mixture,  
Features a screen for making necessary adjustments,  
Wide oscillation range and adjustable mixing time (from 1 to 99 seconds),  
Mixing frequency: 4500 rpm,  
Safe use with automatic shut-off function when the front is open,  
Easy to use, quiet and easy to clean,  
Includes:  
Power cord,  
Instruction manual in Polish, English,  
Dimensions: 35 x 25 x 20 cm,  
Weight: 3 KG,  
Has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

