



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **VDW R-5 RAYPEX 5 endometer (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
German production ,  
Power supply: 230 V,  
Emergency power supply: via 1 built-in 3.6 V / 750 mAh battery,  
Frequency: 50/60 Hz,  
With color visualization and magnification of the last apex section  
RAYPEX 5  
accurately identifies the beginning of the physiologic apex and shows where the anatomic apex is. Illustrates "live" movement of the instrument along the entire root canal - optimal clinical results with a 95% probability of determining the point between the physiological apex and the anatomical apex,  
Measuring the length of the canal with the endometer is completely unobtrusive for the patient and causes no pain. In addition, the measurement can be repeated many times during the procedure, so that the canals can be properly prepared without having to take radiographs several times,  
Features:  
Unique visualization of the apical constriction area,  
Display calibrated to the apical constriction area, separate calibration not necessary,  
"Live" presentation of file movement along the entire canal,  
Resistance measurement in a multi-frequency measurement system,  
Optimized software to achieve even better measurement accuracy and efficiently manage battery power,  
Screen: 5" inch liquid crystal display,  
Included:  
AC adapter with power cord,  
Touch electrode,  
2x attachments,  
Plug tester,  
Device dimensions: 13 x 10.5 x 4 cm,  
Weight: 0.5 KG,  
Has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

