



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DOCTOR'S FRIEND surgical suction machine DF-506B (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

Taiwanese production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Surgical suction machine

DF-506B

is an ideal choice for both hospitals and private health centers. It is a suction unit providing high vacuum and a maximum flow of up to 40 l/min,

DF-506B

is equipped with a modern, maintenance-free and also oil-free pump, thanks to which it provides low noise levels and does not emit harmful oil vapors. The construction made of aluminum alloy guarantees low weight, making the device lightweight as well as comfortable to transport,

Specifications:

Maximum flow rate up to: 40 l/min,

Maximum pressure: 0 ~ 720 mm/Hg,

Suction machine retrofitted with autoclavable jar (121° C) SERRES with 2 l. capacity, including disposable cartridge,

Includes:

Power cord,

Autoclavable jar (121° C.) SERRES, 2 l. capacity, with disposable cartridge,

Overall dimensions: 33 x 30 x 50 cm,

Weight: 12 KG,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

