



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Video-gastroscope PENTAX EG-290Kp (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Japanese production, (2017),
Technical specifications:
Viewing direction: straight,
Viewing angle: 140°,
Depth of field: 5-100 mm,
Working part diameter: 9.8 mm,
Operating channel diameter: 2.8 mm,
Maximum tip bending angle:
Top: 210°,
Bottom: 120°,
Right: 120°,
Left: 120°,
Working length of the sighting part: 1050 mm,
Total length: 1373 mm,
Included:
English instruction manual,
Weight of set: 4 KG,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

