



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDELA DOMINANT 50 surgical suction device (Reconditioned)

Product temporarily unavailable

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

Swiss production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 150 W,

Dominant 50

is a high-quality surgical suction machine that guarantees maximum suction performance and quiet, reliable operation,

Product Features:

Maintenance-free low-speed piston pump,

Use in procedures requiring rapid suction of large volumes of fluids,

Mainly recommended for operating rooms, neurosurgery, plastic surgery and general surgery,

Piston-cylinder technology ensures that the desired vacuum is reached in seconds and that the suction capacity is high,

No vibration due to low rotational speed,

Adapted for continuous operation (24 h / day),

Long life of the device,

Infinitely adjustable suction power,

Built-in pressure gauge,

Diaphragm vacuum regulator allows precise vacuum setting,

Located on the side of the pump, the safety kit prevents overflow and reduces the risk of contamination,

Mobile base has a rail for mounting jars and accessories,

Mounted on the rail with bones,

Possibility to use jars of different capacities,

Specifications:

Capacity: 50 L / min,

Infinitely adjustable suction power: max 90 kPa (675 mm Hg),

Time to obtain vacuum from 0 to 85 kPa on 2 L tank: 14 s,

Included:

User manual in EN, DE, FR, IT, NL, SV, FI, NO, ES, PT (.PDF)

2 x autoclavable jar (121° C) SERRES 2 L. with disposable cartridge,,

Suction hose with a length of 150 cm,

Power cord,

Dimensions:

Suction jar: 39 x 25 x 33 cm,

General: 48 x 48 x 108 cm,

Weight: 24 kg,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for a period of 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

