



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PANASONIC WV-CM1780 medical monitor (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Philippine production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 73 W,
Specifications:
Display type: CRT,
Screen diagonal: 17" inches,
Inputs: S-VIDEO, AUDIO RCA 1, AUDIO RCA 2, RJ45,
Outputs: AUDIO RCA 1, AUDIO RCA 2, RJ45,
Selectable menus in ENG, FR, GER, ESP,
Included:
Power cord,
Dimensions: 40.5 x 35 x 42.5 cm,
Weight: 18.5 KG,
Has a current inspection and is ready to work,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

