



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Surgical suction machine NEW HOSPIVAC 350 CA-MI (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Italian production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 230 VA,  
NEW HOSPIVAC 350 suction machine  
is used to remove various secretions from internal organs. Its great advantage is that it works continuously and does not lose efficiency.  
Thanks to the fact that it works quietly, it can be used in the operating room or ICU,  
The device is enriched with a special pump that assists in suction. Thanks to the locking mechanism on two wheels, the product stands stably in place,  
Technical specifications:  
Maximum vacuum: -90 kPa,  
Maximum flow rate: 60 l./min,  
Noise level: 51.7 dBA,  
Duty cycle: continuous operation,  
Included:  
Power cord,  
2x 2-liter jar,  
Wired foot switch,  
Dimensions: 46 x 42 x 85 cm,  
Weight: 13 KG,  
It has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

