



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **HUNTLEIGH AKRON diagnostic and surgical chair (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
U.K. production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Upholstery with medical approval in blue color (it is possible to replace the upholstery with a new one, from our color pattern for an additional fee),  
HUNTLEIGH AKRON chair - is designed for treatment and diagnostic examinations for a wide range of disciplines, from simple examinations to heavier manipulative and manual therapeutic techniques,  
3 - segment recliner,  
Steel construction in white color,  
Equipped with:  
Adjustable armrests on both sides,  
Divided leg segment,  
Paper roll holder with a width of 46 cm,  
Electric height adjustment of the table by foot attachment in the range: from 43 to 93 cm,  
Manual adjustment of:  
Backrest angle in the range of 74°,  
Angle of the footrests within the range: 34° / 0° / -90°,  
Mobile chair with four casters with built-in braking feet,  
Dimensions:  
Backrest: 64 x 82 cm,  
Seat: 64 x 52 cm,  
2 - segment leg rest: 58 x 37 cm,  
Total length in reclining position: 180 cm,  
Safe working load: 225 KG,  
Maximum lifting capacity: 190 KG,  
Weight: 90 KG,  
Has a current inspection and is ready for use,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

