



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



HUNTLEIGH AKRON rehabilitation chair (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
U.K. production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Upholstery with medical approval in blue color (it is possible to replace the upholstery with a new one, from our color pattern for an additional fee),
HUNTLEIGH AKRON chair - is designed for treatment and diagnostic examinations for a wide range of disciplines, from simple examinations to heavier manipulative and manual therapeutic techniques,
3 - segment recliner,
Steel construction in white color,
Equipped with:
Adjustable armrests on both sides,
Divided leg segment,
Paper roll holder with a width of 46 cm,
Electric height adjustment of the table by foot attachment in the range: from 43 to 93 cm,
Manual adjustment of:
Backrest angle in the range of 74°,
Angle of the footrests within the range: 34° / 0° / -90°,
Mobile chair with four casters with built-in braking feet,
Dimensions:
Backrest: 64 x 82 cm,
Seat: 64 x 52 cm,
2 - segment leg rest: 58 x 37 cm,
Total length in reclining position: 180 cm,
Safe working load: 225 KG,
Maximum lifting capacity: 190 KG,
Weight: 90 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

