



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Gynecological chair FAMED ŻYWIEC FG-06.0 (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- Polish production,
- **Medically certified upholstery in green (the upholstery can be replaced with a new one from our colour chart for an additional fee),**
- The structure of the armchair is made of powder-coated stainless steel in white,
- **Technical specifications:**
 - Manual backrest adjustment in the range of: 70°,
 - Manual adjustment of seat tilt in the range: 22° / 0° / -10°,
 - Fixed seat height: 82 cm,
 - Possibility of setting the seat in the Trendelenburg position,
- **It has:**
 - Double-sided rails for mounting additional equipment,
 - Manually adjustable bilateral knee rests,
 - Patient handles attached on both sides,
 - Pull-out waste bowl measuring 30 x 33 x 11 cm,
 - Dedicated patient level,
 - Paper hanger 40 cm wide,
- **Included in the set:**
 - Patient degree,
 - 2x patient knee rest,
 - 2x patient holder,
 - Pull-out waste bowl,
- **Dimensions:**
 - Backrest: 56 x 85 cm,
 - Seat: 56 x 38 cm,
 - Total length in lying position: 125 cm,
 - Total: 90 x 110 x 159 cm,
- Weight: 80 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for domestic market (Poland),**
 - **3 months to international market,**
 - **Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,**
- Financing options (Poland only): Installments, Leasing, Loan,

