



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Gynecological chair SCHMITZ MEDI-MATIC 115.305 (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- Power supply: 230V,
- Frequency: 50/60 Hz,
- Power: 400 W,
- Medically certified upholstery in dark grey (**the upholstery can be replaced with a new one from our colour chart for an additional fee**),
- **Schmitz Medi - Matic** is used for gynecological and urological examinations and procedures,
- **Technical specifications:**
 - Possibility of electromechanical adjustment using a foot switch:
 - Seat height range 84.5 - 106 cm,
 - Synchronized adjustment of the backrest and seat,
 - "Zero" position function,
 - Smooth manual adjustment of the backrest by hand, using a gas spring in the range: from 0° to 45°,
 - Double-sided Goepel knee pads adjustable in many dimensions and attached to the railing,
 - 43 cm wide paper mount,
 - Backrest with detachable middle segment,
- **Included in the set:**
 - Power cord,
 - Wired foot controller,
 - Goepel knee pads adjustable in many dimensions and mounted on side rails – 2 pcs.,
 - Retractable patient step,
 - Pull-out waste bowl with dimensions: 32 x 30 x 6 cm,
 - Double-sided patient handles,
- **Dimensions:**
 - Seat: 59 x 31.5 x 12.5 cm,
 - Backrest: 59.5 x 87 x 12.5 cm,
 - Footrest: 43.5 x 33.5 cm,
 - Knee rests: 28.5 x 18.5 cm,
 - Base: 53.5 x 79.5 x 5 cm,
 - Overall dimensions: 87 x 128.5 x 87/168 cm,
- Maximum load capacity: 135 KG,
- Weight: 161 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for domestic market (Poland),**
 - **3 months to international market,**
 - **Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,**
- Financing options (Poland only): Installments, Leasing, Loan,

