



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **Gynecological chair SCHMITZ MEDI-MATIC 115.755 (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- Frequency: 50/60 Hz,
- **The SCHMITZ MEDI-MATIC** gynecological chair is used to perform gynecological and urological examinations and procedures,
- **Medically certified upholstery in orange (the upholstery can be replaced with a new one from our colour chart for an additional fee),**
- **Technical specifications:**
  - **Smooth electric adjustment via foot controller:**
    - Heights range: 60 - 90 cm,
    - Backrest adjustment angle in the range:  $\pm 52^\circ$ ,
    - Seat up/down angle:  $\pm 5^\circ$ ,
    - Trendelenburg position:  $19^\circ$ ,
  - **Smooth electric adjustment via handheld remote control:**
    - Heights range: 60 - 90 cm,
    - Backrest adjustment angle in the range:  $\pm 52^\circ$ ,
    - Seat up/down angle:  $\pm 5^\circ$ ,
    - Trendelenburg position:  $19^\circ$ ,
    - Double-sided footrests in the range of:  $\pm 45^\circ$ ,
- Quick change of seat position using the foot controller,
- The low starting position of the seat allows easy and safe entry without the need to use a step,
- Paper roll storage,
- 48 cm wide paper roll hanger,
- Pull-out waste bowl from under the seat section,
- Gentle start and ergonomic movement when changing positions,
- **Included in the set:**
  - Power cord,
  - Wired foot controller,
  - Wired remote control,
  - A bowl for waste,
  - Headrest,
- **Dimensions:**
  - Bases: 105 x 58 cm,
  - Waste bowls: 31 x 26 cm,
  - Overall: 150 x 135 x 86 cm,
- Maximum load capacity: 200 KG,
- Weight: 165 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
  - **6 months for domestic market (Poland),**
  - **3 months to international market,**
  - **Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,**
- Financing options (Poland only): Installments, Leasing, Loan,

