



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **MIDMARK PROMOTAL 3181 dialysis chair (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- French production,
- Power supply: 230V,
- Frequency: 50/60 Hz,
- Power: 230 W,
- **Medically certified upholstery in grey (the upholstery can be replaced with a new one from our colour chart for an additional fee),**
- **The MIDMARK PROMOTAL 3181 dialz chair provides free adjustment with the possibility of positioning in a wide range,**
- The chair has excellent mobility, making it suitable for free movement around medical facilities,
- **Technical specifications:**
  - **Adjustment via wired handheld remote control:**
    - Seat height adjustment range: 56 – 80 cm,
    - Backrest adjustment range: -7° / 0° / 74°,
    - Leg rest adjustment range: -57° / 0° / 5°,
    - Tilt to Trendelenburg position,
    - The remote control has a programmed button to set the seat to its initial position,
  - **Mobile chair on 4 wheels with built-in brake,**
  - **Bilateral freely adjustable patient armrests,**
  - **Freely adjustable drip hanger with the possibility of double-sided mounting,**
- **Included in the set:**
  - Patient headrest,
  - 2 pcs. - armrest,
  - IV drip hanger ,
- Dimensions:
  - Backrests: 58 x 80 cm,
  - Seats: 58 x 53 cm,
  - Leg rests: 58 x 46 cm,
  - Total: 95 x 115 x 145/169 cm,
- Maximum load capacity: 135 KG,
- Weight: 95 KG
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
  - **6 months for domestic market (Poland),**
  - **3 months to international market,**
  - **Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,**
- Financing options (Poland only): Installments, Leasing, Loan,

